

Step 1 – Order Header & Audiogram

- Please fill in client / Header details
- Enter delivery date required
- Indicate the experience level
- Enter Audiometric data (mandatory for AOV)
- Select left/right side

Step 2 – Performance level

- Select performance level :
B90 Premium, B70 Advanced
- Select the user controls

Step 3 – Product Options

- Select the custom product options
- Use the track «smallest possible» for fast Ordering
— *Smallest possible*
- The shaded box options will increase device size

AOV = Acoustically Optimized Venting –
Phonak will set the correct vent diameter and style in relation to Audiogram / shell geometry / feedback canceller and required low frequency gain – this is the recommended vent type.

AOV-O = Acoustically Optimized Venting - Open
This specific venting is more open and recommended for **first time users** and clients who are sensitive to "occlusion".

- Please let us know about special wishes
- Select warranty and rush order details

Cosmetic aspect

Human ear canals are not fully symmetrical. Therefore binaural custom in-the-ear devices cannot always be build in the same size and position in the ear.

As an option we can design the two devices to look similar in the ear. As a consequence one device may not be as small as possible.


Please request this with the «symmetrical build» (SY) check box.

Shell Styles




Power Level

Virto B-Titanium is designed for smallest device size. Powerlevel P/SP does increase devices size AOV-O to be used for hearing loss below 40dB only (see graph)



Impression taking



Recommendation for taking impressions:

- Use EasyView Otoblock
- Use Dreve Otoform A softX 25 shore A (#2956)
- Always fill the complete concha with impression material

